


Monthly alarm health check

Performing regular tests of the LiveLife Mobile Alarm

Like other electronic devices, it is important to conduct regular checks of an alarm to make sure it is working correctly. We suggest a monthly health check which includes the following:

- 1 Check that the **"green"** light on your alarm is flashing once quickly every 3 seconds. This means you have a connection to the mobile network.



- 2 Ask a relative or friend to **call** your alarm and you can answer in hands-free, two-way voice.



- 3 Hold down the **SOS button** for one second until you hear the voice prompt. Click the middle SOS button again to cancel the emergency call and text sequence.



- 4 Check your alarm for any obvious signs of damage. Make sure the three brass discs on the back of the alarm are bright and shiny.



Lifetime support

Your alarm comes with lifetime support. If you have any concerns, give us a call on 1800 936 774 and press 2 for Customer Support. 8.30am to 5.00pm, Mon-Fri, AEST.

Support videos

A video that demonstrates how to do a family test run of your mobile alarm can be viewed at:

www.livelifealarms.com.au/videos